arzum

SATEEN TOUCH

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Üretici:

Foshan Jindi Electric Appliance Co., Ltd

No. 13 Baiye Avenue, Xiqiao Science and Technology Industrial Park, Xiqiao
Town, Nanhai District, Foshan City, Guangdong Province, China
0086 75786868180

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
0850 222 1800

www.arzum.com

AR5178-02 ARZUM SATEEN TOUCH RECHARGEABLE IPL HAIR REMOVAL SYSTEM

- 1- UV-filtered light window
- 2- Power/Mode button
- 3- Energy level increase button
- 4- Energy level decrease button
- 5- LED display
- 6- Cooling function indicator
- 7- Energy level indicator
- 8- Charging indicator
- 9- Automatic mode indicator
- 10 - Body / Bikini / Face mode indicators

- 11- Light pulse button
- 12- Air vents
- 13- Magnetic face attachment
- 14- Magnetic bikini attachment
- 15- Skin contact sensor and indicator
- 16- Adapter port
- 17- Adapter
- 18- Cleaning cloth
- 19- Travel pouch

Thank you for purchasing **Arzum IPL Hair Removal System**. This is one of the products designed to make life easier for you.

Before using the device, carefully read the instructions for use. Keep the user manual, warranty certificate, receipt, or invoice, and if possible, the box as well. If you give this device to someone else, please give them the user manual along with it.

SAFETY INSTRUCTIONS

- Use the device only for removing unwanted hair. After prolonged use, a reduction in the number of growing hairs may be observed. Use for any other purpose is inappropriate and may be dangerous.
- Do not use on animals.
- Do not use the device on any surface other than the skin.
- Do not use the device for commercial purposes.
- When not in use, during cleaning, or in case of malfunction, turn off the device and unplug it.
- If you will not be using the device again, it is recommended that you render it unusable by removing all dangerous parts, especially for the safety of children who may be tempted to play with the product.
- Provided that training is given on safe use and potential hazards, the device can be used by children aged 8 and above, as well as individuals with physical, sensory, or mental disabilities, or those lacking experience and knowledge. However, children should not play with the device; they should not be allowed to clean or maintain the device on their own.
- If your device is damaged for any reason, do not use it, do not open it, and do not attempt to repair it yourself. In case of any malfunction, take

it to an **Arzum Authorized Service Center**.

- For safe and efficient operation of the device, use only original parts and adapters.
- Do not use products with damaged cables. Contact an Arzum Authorized Service Center.
- Do not allow children to play with the product packaging. Plastic bags should not be left where children can reach them due to the risk of suffocation.
- Before connecting the device, ensure that the main voltage information matches the voltage information specified in the product's technical specifications. The technical specifications are listed on the device and on the adapter.
- Protect the device from sunlight, rain, and moisture.
- Do not use the device near water or on wet surfaces (bathtub, shower, sink, etc.).
- Do not use the device if it gets wet. Using a wet device may cause electric shock.
- Do not use the device if the light window (1) is cracked or missing.
- Do not touch the light window during or immediately after use; it may be hot.
- Do not use petroleum-based or flammable cleaning agents (alcohol, acetone, etc.); they pose a fire hazard.
- Do not cover the cooling vents. The device may overheat, posing a risk of electric shock or fire.

General health status and chronic illnesses

Do not use the epilator device if you have any of the following health conditions, or consult your doctor beforehand:

- If you are under 18 or over 65 years of age.
- If you have diabetes or any other systemic or metabolic disease.
- If you have congestive heart disease.
- If you have active implants such as a pacemaker, neurostimulator, or insulin pump.
- If you have had a disease that suppresses the immune system (e.g., HIV infection or AIDS).
- If you are taking medications that suppress the immune system.

- If you have any clotting disorder.
- If you have photosensitive epilepsy.

Sensitivity and Skin Sensitivity Conditions

- If you have a skin tone that is not suitable according to the **"Suitable skin tones"** table in the user manual.
- If you have a condition related to light sensitivity, such as polymorphic light eruption (PMLE), solar urticaria, or porphyria.
- If your skin is sensitive to light and easily develops rashes or allergic reactions.
- If you are using medications or substances that increase light sensitivity and their package inserts indicate a risk of photoallergic or phototoxic reactions.
- If you are using painkillers that reduce the skin's sensitivity to heat.

Medical History

- If you are pregnant or breastfeeding.
- If you have used Accutane, Roaccutane, or any other form of isotretinoin within the last 6 months.
- If you are undergoing skin treatment with Alpha Hydroxy Acid (AHA), Beta Hydroxy Acid (BHA), topical tretinoin, or azelaic acid, or if such treatment was administered within the last 7 days.
- If you are taking blood-thinning medications, including high doses of aspirin, and require at least a 1-week washout period for these medications.
- If you have undergone radiation therapy or chemotherapy within the last 3 months.

Skin Conditions and History

- If you have a history of collagen disorders, such as keloid formation or slow wound healing.
- If you have active skin cancer, a history of skin cancer, or a history of local cancer in the areas to be treated.
- If you have a skin pigmentation disorder, such as vitiligo.

Conditions Related to the Treatment Area

- If you have a history of vascular disease in the treatment area, such as varicose veins or vascular dilation (vascular ectasia).
- If there is an infection, eczema, burn, folliculitis, open wound, abrasion,

cold sore (herpes simplex), lesion, or bruise (hematoma) in the areas where the application will be performed.

- If the skin is irritated (redness, cuts), sunburned, tanned, or if a tanning product has been applied.
- On moles, freckles, large veins, dark pigmented areas, scars, or skin imperfections without consulting your doctor.
- On areas with warts, tattoos, or permanent makeup.
- Around the eyes and eyebrows. It can only be used on the body and on the female face (under the cheekbones).
- On the lips, scalp, nipples, between the legs, around the nipples, labia minora, vagina, anus, nostrils, and inside the ears.
- Areas where long-lasting deodorant is used. This may cause skin reactions.
- On or near artificial objects such as silicone implants, subcutaneous injection channels (e.g., insulin dispensers), or piercings.

For men: Shoulders and below (chest, back, arms, abdomen, and legs)

WARNING: Do not use on the face, neck, or genital area. Use on facial hair or beard may cause permanent or uneven results; this may lead to changes in physical characteristics.

NOTE: This list does not cover all possibilities. If you are unsure whether you can use the device, be sure to consult your doctor.

IS IPL HAIR REMOVAL TECHNOLOGY SUITABLE FOR YOU?

For IPL (Intense Pulsed Light) technology to be effective, there must be a clear contrast between skin tone and hair color. The best results are achieved by users with fair skin and dark hair. IPL light is absorbed by melanin, the natural pigment found in skin and hair. This absorption heats the hair follicle, rendering it ineffective.

However, effectiveness may decrease or risks may arise in the following cases:

Red, white, gray, and light yellow hair may not absorb enough light due to their low melanin content, and the desired results may not be achieved.

Very dark skin tones may absorb too much light due to high melanin content, which can cause burns on the skin.

NOTE: For first-time users, energy level 1 is recommended to assess skin sensitivity. After completing the skin test, wait for 24 hours and

check for any reactions. Higher energy levels provide more effective results. For subsequent applications, it is recommended to use the device at an increased energy level.

SKIN TONE & HAIR COLOR CHART TREATMENT SCHEDULE

Please check your skin tone and hair color against the table before using the device to ensure compatibility. Using the device with incompatible combinations may cause skin irritation, burns, or ineffective results.

SKIN COLOR						
HAIR COLOR						

APPLICATION SCHEDULE

Recommended frequency of use for faster and more effective results:

Weeks 1–2: Apply 3 times per week.

Weeks 3–4: Continue applying 2 times per week.

Weeks 5–8: Apply 1–2 times per month.

NOTE: After the 8th week, if additional hair growth is observed, maintenance is recommended every 1–3 months.

The interval between applications may vary depending on the hair growth rate and the body area; effectiveness may differ from person to person.

Visible Results After 8 Weeks:

Week 1		You will still see some hairs growing. These are likely to be hairs that were missed during treatment.
Week 2		You should see a reduction in hair growth. However, many hairs may still not have been treated in their growth phase.
Week 4		You should see a significant reduction in hairs within the treated area.
Week 8		Any remaining hairs should be finer and lighter in color.

Cooling Function

The cooling technology in the fixed head cools the skin, making the application more comfortable. It is automatically activated when the device is turned on, and the (❄️) icon appears on the screen. The face and bikini heads do not have a cooling function.

Skin Tone Analysis Sensor

Automatically adjusts the pulse intensity by detecting skin type and tone. This allows the device to select the most appropriate light level, ensuring a safe and effective application.

The skin tone analysis sensor works only with the fixed body head. Before attaching the head, detect the skin tone of the area you want to treat. Then, attach the bikini or face head and perform the application according to the detected level.

Skin Contact Sensor

Detects whether the device is properly positioned on the skin. When the sensor is activated, the rear LED indicator lights up white, indicating the device is ready to emit a light pulse. This helps prevent unnecessary flashes.

Magnetic Heads

To achieve the best results, it is important to use the correct head designed for the specific body area you want to treat.

In addition, thanks to their magnetic structure, the heads can be easily attached and removed, providing practicality and comfort during use.

This product includes two magnetic heads:

Bikini Area Head

· Provides effective and safe application on the bikini area, which has a more sensitive and dense hair structure. It offers comfortable use by establishing more controlled contact with the skin.

Face Area Head

· Specially designed for sensitive areas such as the face. Its small application surface allows for precise and safe light emission on areas like the upper lip, chin, and cheeks.

· Since these heads are used on different body areas, it is important to use them separately for hygiene reasons.

USE OF THE DEVICE

Charging the Device

· Before using the device for the first time, fully charge the battery if it is completely empty.

· It takes approximately 2 hours for the battery to fully charge. A fully charged battery provides at least 260 flashes at light intensity setting 5.

· Plug the adapter into a power outlet.

· During charging, the battery icon will blink; once fully charged, the icon will remain steadily lit.

· After charging is complete, unplug the adapter from the outlet.

· When the battery level is low, the device automatically turns off. In this case, you can recharge the device or continue using it in wired mode.


CAUTION: Do not cover the adapter or the device during charging.

NOTE: It is normal for the adapter and the device to heat up during charging.

This device is equipped with overheating protection; if the ambient temperature exceeds 70 °C, the charging process will not start.

Preparation Before Application

- Prepare the skin by shaving the hair in the area to be treated, then clean and dry the area.
- This step allows the IPL light to reach the hair follicles directly and increases the effectiveness of the application. Methods other than shaving may remove the hair follicle, which can negatively affect the performance of the application. If hair remains on the surface, a burning smell will occur during the application. Therefore, it is not recommended to remove hair with gel, cream, or an epilator.

Operating the device

- The device can be used both wired and wireless.

For wired use: Connect the IPL device to the adapter and plug the adapter into the power outlet. Press the power button to turn on the device. This mode provides uninterrupted power for long sessions.

For wireless use: Press the power button to turn on the device. This mode offers freedom of movement and practical use.

Operating Modes

- The device has three different modes for different needs: Application, Manual, and Automatic.

Application mode: One of three different modes with pre-set light intensity and pulse rate according to the area to be treated can be selected. These modes ensure that the device works appropriately for the relevant area.

Manual Mode: In this mode, the device does not emit light pulses automatically; you must press the flash button for each light pulse. It is suitable for small and sensitive areas such as fingers and the upper lip.

Automatic Mode: In this mode, the device emits light pulses automatically. Press and hold the flash button for 2 seconds to activate or deactivate the automatic mode. The () icon will appear on the display.

This mode is ideal for large treatment areas such as arms and legs.

Skin Analysis Sensor: To activate the skin analysis sensor, set the energy level to "0". The device will automatically adjust the energy level according to the skin tone of the treatment area.

- Press and hold the power button for 2 seconds to turn on the device. The fan will start operating automatically.

- Your device has 5 different energy levels. After selecting the desired level, you can perform manual flashes.
- Place the light output window (1) firmly against the skin on the area to be treated. When the window is in full contact with the skin, the LED light on the back of the device begins to flash and stays on steadily when full contact is achieved. This indicates that the device is ready to emit light. Press the button to emit a light pulse. You may feel a slight warmth on the skin.

Move the light window to the next area and repeat the process.

- To use the automatic mode, press and hold the flash button for 2 seconds.
- To start treatment, place the device vertically against the skin. Then glide the light window across the skin; the device will automatically emit light pulses when it detects contact.

In this way, there is no need to press the flash button repeatedly.

To deactivate the function, press and hold the flash button for 2 seconds again; the () icon on the screen will disappear.

- Make sure that the entire area has been treated, but do not flash the same spot twice (avoid overlapping light pulses).


NOTE: It is normal for the device to become slightly warm during use.

CAUTION: If the device is not used for 3 minutes, it will automatically shut down.

If the internal temperature reaches 75°C, an error message will appear on the screen.


Please allow the device to cool down in a ventilated area; it will resume normal operation afterward.

- When the treatment is complete, press and hold the power button to turn off the device.
- The fan will continue to run for approximately 5 seconds to cool the device. During this time, a countdown timer will appear on the screen, and the device will automatically power off afterward.

CLEANING

- Before cleaning the device, make sure it is turned off and unplugged.
- The light window may become hot during use; ensure it has cooled completely before cleaning.
- Do not use dish sponges, abrasive cleaning materials, gasoline, acetone, or other liquids that may damage the surface.
- After use, wipe the light window with the special cleaning cloth provided with the product.
- If necessary, wipe the exterior of the device with a clean, soft cloth.
- To get the best performance from the device, clean it after each use.

SIMPLE TROUBLE-SHOOTING

When fault occurs during the usage of the instrument, users can follow the following guidance to eliminate some simple faults. If the problem can not be solved after all steps are finished, please contact the Customer Service Department so as to solve faults as soon as possible.

Problem	Solution
The device does not power on when the Power button is pressed.	· Check the device's battery level. · If the battery is completely drained, charge the device and try again. · If you are using the device with a cable, make sure the adapter is properly connected to both the device and the power outlet.
The device does not emit a light pulse when the Light Pulse button is pressed.	· Check whether the device has overheated. The device will not emit pulses if it is too hot. · Ensure that the treatment window is in full contact with the skin and that the Light Pulse READY indicators on both sides of the device are illuminated.
The energy level does not change when the Energy Level button is pressed.	Pressing the Energy Level button repeatedly or rapidly may trigger the device's 3-minute safety lock. Wait for 3 minutes, then press the button again to select the desired energy level.
Significant discomfort is experienced on the skin during device use.	• Hair trimming: Shave the treatment area prior to use to shorten hair. • Recent tanning: Exposure to sunlight or artificial tanning products shortly before use may increase skin sensitivity. Avoid sun exposure and tanning devices for at least 2 weeks prior to treatment. • Energy level adjustment: Reduce the energy level to a comfortable setting.
The device is powered on, lights are illuminated, but no pulse is emitted.	The device's safety system may have activated due to overheating. Allow the device to cool; normal pulsing will resume. This is a safety feature and does not indicate a malfunction.

Model No	AR5178
Power	48W
Light Output Area (Spot Size)	3.3 cm ²
Wavelength	600-1200 nm
Total Energy	20J
Energy Density	6.06 cm ²
Repetition Rate	1.1-4 Seconds (depending on the selected energy level)
Operating Mode	Auto Mode / Manual Mode
Product Size	15.34*8.17*19.05 cm
Additional Header (Spot Size)	1 cm ² / 2 cm ²
Battery	7.4Vd.c.; 1500 mAh; 11.1Wh

Product Life 7 years

Technical Specifications

Main Unit: 24V $\equiv$ 2A

Adapter Model: DZ048CHL240200V

Adapter Input: 100-240V~50/60Hz 1.5A

Adapter Output: 24V $\equiv$ 2A 48.0W

HANDLING AND SHIPMENT

- During handling and shipment, your appliance must be kept in original package in order to prevent damages to its parts.
 - Keep at normal position during shipment.
- Do not drop the product when transporting and protect from impacts.
- After delivery of the product, failures and damages that may occur while transporting are not covered by warranty.

The product is in compliance with the WEEE Directives.

WEEE

This product contains recyclable materials that compliance with the WEEE directives. Do not dispose this product as unsorted municipal waste. Please contact your local municipality for the nearest collection point.

Manufacturer:

Foshan Jindi Electric Appliance Co., Ltd

No. 13 Baiye Avenue, Xiqiao Science and Technology Industrial Park, Xiqiao Town, Nanhai District, Foshan City, Guangdong Province, China
0086 75786868180

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No. 13 Baiye Avenue, Xiqiao Science and Technology Industrial Park, Xiqiao Town, Nanhai District, Foshan City, Guangdong Province, China

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
0850 222 1 800

www.arzum.com

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are approximately 20 lines visible. The top-left corner of the paper is rounded. The paper appears to be a standard notebook page.

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