arzum

SENSIBLE TOUCH

AR5300-02 ŞARJLI AĞIZ DUŞU
AR5300-02 RECHARGEABLE WATER FLOSSER

KULLANMA KILAVUZU
INSTRUCTION MANUAL
7.4 W

2 YIL YEARS STANDART WARRANTY

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
AR5300-02
ARZUM SENSIBLE TOUCH
RECHARGEABLE WATER FLOSSER

1- Spray Nozzle
2- Nozzle Holder
3- Nozzle Storage Cover
4- On/off swich
5- Mode Button
6- Mode Indicator
7- Charging Indicator
8- Reservoir
9- Reservoir Lid and Mirror
10- USB Charging Cable
11- Charging Port

Thank you for purchasing **Arzum** brand flosser. We wish this product to make life easier for you, as with other **Arzum** products offered to you.

In order to obtain utmost output from your appliance, please read carefully this user's manual and save it for further use.

IMPORTANT SAFEGUARDS

· Make sure that the voltage indicated on **Arzum Sensible Touch Flosser** corresponds with the main voltage in your home. Our company shall not be liable for any failures owing to any discrepancy and such failures shall not be covered within the guarantee.

· Your device is for household use only, not for commercial or industrial usage. Otherwise, inconvenient uses shall not be covered by the guarentee.

· This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensor your mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance.

· Cleaning and user maintenance shall not be made by children without supervision.

WARNING: Keep the appliance dry.

The appliance is only to be used with the USB cable provided with the appliance.

This appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance.

· Consult your doctor before using this device to obtain medical information.

· If you have had dental or gum surgery in the past 2 months, consult your dentist before using this device.

· If you use a pacemaker or any other surgically implanted device, consult your doctor before using the oral irrigator.

· This product is designed only for cleaning teeth and gums. Do not use it for other purposes (e.g., eyes, nose, or ears).

· If bleeding occurs after use, if bleeding continues after one week of use, or if you feel discomfort or pain, stop using the product and consult your dentist or doctor.

· This device is a personal care product and should not be used on multiple patients in dental clinics or hospitals.

· Do not use the product without the nozzle attached.

· Do not apply excessive pressure to the tip of the nozzle.

· Replace the nozzle every 6 months or earlier if signs of wear appear.

· Do not use the device with accessories other than those specifically designed for it.

· Do not attempt to disassemble any part of the device or insert metal objects into it.

· Do not disassemble or modify the device.

· Use the device only as recommended and for its intended purpose.

· Do not leave the device unattended while it is plugged in.

· Do not use the device if it has been dropped or if there are visible damages.

· In case of malfunction, take the device to an authorized **Arzum service** center.

· Always use only original spare parts with your device.

· Keep the charger away from water.

· Do not fill the water tank with hot water. Hot water can cause burns.

· To prevent physical injury, do not fill the tank with hot water.

· Do not put foreign objects into the water tank.

· Do not clean any parts of the product in a dishwasher.

· Do not place the product or its batteries in a microwave oven or on induction cooktops.

· Charge and use the device within a temperature range of 5°C to 40°C (41°F to 104°F).

· Do not charge the product outdoors or near hot surfaces.

· Keep the product and its batteries away from fire, direct sunlight, and high temperatures.

· If the product overheats, emits odor, changes color, or takes longer than usual to charge, stop using it and contact **Arzum authorized service.**


· To prevent the battery from overheating or releasing

toxic or hazardous substances, do not puncture, disassemble, modify, or damage the product or the battery.

· Do not short-circuit, overload, or reverse-charge the batteries.

· Avoid contact of damaged or leaking batteries with skin or eyes. If contact occurs, wash the affected area thoroughly with water and consult a doctor immediately.

CHARGING THE DEVICE

· When charging the product for the first time, make sure the battery is fully charged. The standard charging time is ≤ 2.5 hours.

· Carefully open the charging port cover (11) located at the top of the device. Be careful not to force the cover backward.

CAUTION: The package includes only a USB charging cable. It is recommended to use a phone adapter with a 5V-1A output to charge the device.

· Connect one end of the charging cable to the charging port on the device, and the other end to a power source with a USB output (such as an adapter, computer, or power bank).

· If the red indicator light is flashing → This means the device needs to be charged.


· If the red light stays on continuously during charging → This indicates the device is currently charging.

· When the green light stays on continuously → This means the device is fully charged.

· After the device is fully charged, close the charging port cover to prevent water from entering and to maintain the device's performance.

USING THE DEVICE

PREPARING FOR USE


· Hold the water tank steady and pull it upward from the top; you will hear a distinct "**click**" sound when the tank is properly positioned. It is important to hear this sound—otherwise, the tank may not be fully secured and could shift back, causing water to leak.

· Open the tank lid by pulling it to the left. Fill the tank with clean water.


· Carefully close the tank lid. Make sure the lid is properly seated when closing.

· When the lid is securely closed, you will hear a "**click**" sound. This indicates that the lid is safely locked in place.


CAUTION: Consult your doctor before using mouthwash or similar solutions. This product is compatible with such solutions; however, the device must be thoroughly rinsed after each use. Inadequate rinsing may cause clogging.

· Remove the nozzle from the storage compartment by pulling it out. Insert the nozzle into the device and turn it clockwise. When you hear a "**click**" sound, it means the nozzle is securely attached. To remove it, turn the nozzle counterclockwise and place it back into the storage compartment.

Power On/Off and Mode Switching

Power Button: Press the button briefly to turn the device on or off. (1)

Mode Button: You can change the mode by pressing the mode button briefly, either before turning the device on or during use. (2)

Soft Mode: Suitable for people with sensitive gums, beginners, or those who prefer a gentle cleaning. It sprays water at a lower pressure. (3)

Standard Mode: Designed for general oral hygiene, ideal for cleaning between teeth and removing plaque. It sprays water at medium pressure. (4)

Personal (DIY) Mode: The DIY mode allows you to adjust the water pressure to your personal preference and saves it for future use, offering a comfortable and customized cleaning experience. This setting is configured while the device is off. (5)

· With the device turned off, press the mode button briefly → DIY mode is activated, and the light stays on steadily.

· Press and hold the mode button → The device starts, and the water pressure increases cyclically from low to high.

· Release the button when your desired pressure level is reached → The device stops, and the light remains steady → The setting is now saved. DIY mode remembers your selected pressure level, so you don't need to set it again each time.

Pulse Mode: Provides pulsating water flow for effective interdental cleaning and gum massage.


· After selecting your desired mode, press and hold the Power button for more than 2 seconds to activate Pulse Mode.
· Release the button to stop the mode.

· To restart, press and hold the button again.

STARTING USE

· Before turning on the device, hold the oral irrigator and lean over the sink.

· Place the nozzle inside your mouth and align it perpendicularly to your gums or between your teeth.

· To prevent splashing, keep your lips slightly closed.

· While cleaning, slightly open your mouth to allow water to flow out comfortably.

· For the best cleaning results, it is recommended to start from the molars and move toward the front teeth.

· When inserting the nozzle into your mouth, make sure the water flows at a 90° angle to the gum line.

· Glide the nozzle along the gum line, pausing briefly between each tooth.

· Continue the cleaning process until both the inner and outer surfaces of your upper and lower teeth are completely cleaned.

CAUTION: Do not press the nozzle hard against your gums or teeth. Doing so may cause gum bleeding.

CLEANING AND MAINTENANCE

· After each daily use, first make sure the device is turned off. Then, empty any remaining water from the tank. Turn the device on and let it run for a few seconds without water to remove any residual liquid. This helps prevent bacterial growth and unpleasant odors.


· Minerals in tap water may affect the performance of the device. Once to three times a month, clean the device by adding warm water and one tablespoon of white vinegar to the tank and running it. Repeat the process until the vinegar smell disappears.

· Rinse the nozzle with water and wipe it clean with a cloth. It is recommended to replace the nozzle every 6 months. If the tip becomes deformed, replace it immediately to prevent reduced performance or potential device malfunction. **(A)**

· Clean the main body with a soft, dry cloth. For stubborn stains, use a diluted neutral detergent, then dry the surface with a cloth. **(B)**

· The water tank is designed to be removable for easy cleaning. It is recommended to clean it once a week using diluted neutral detergent and a small brush for detailed scrubbing. **(C)**

· Store the product in a dry, well-ventilated place, out of reach of children.

TROUBLESHOOTING

The problem	Why?	Solution
Power indicator does not light up during charging.	Improper power connection.	Reconnect the power supply correctly.
Short operating time after charging.	· Low ambient temperature. · Insufficient charging time.	· Charge and use the product in an indoor environment between 15°C–35°C. · Charge the product properly.
Device does not turn on.	Battery power is too low.	Charge the product.
Low water pressure.	· Nozzle is dirty or clogged. · Nozzle is damaged.	Replace the nozzle.
Water leakage.	Tank, hose, or nozzle is cracked.	Replace the tank, hose, or nozzle and contact **Arzum Authorized Service.**

Product Life: 5 years

TECHNICAL SPECIFICATIONS

Voltage: 3.7 V DC Power: 7.4 W

Input: Power 5 V- 1A

Battery Capacity 920mAh

IPX7

HANDLING AND SHIPMENT

- During handling and shipment, your appliance must be kept in original package in order to prevent damages to its parts.
- Keep at normal position during shipment.
- Do not drop the product when transporting and protect from impacts.
- After delivery of the product, failures and damages that may occur while transporting are not covered by warranty.

The product is in compliance with the WEEE Directives.

WEEE

This product contains recyclable materials that compliance with the WEEE directives. Do not dispose this product as unsorted municipal waste. Please contact your local municipality for the nearest collection point.

CE

Manufacturer:

NINGBO HAOJIA ELECTRICAL APPLIANCE CO., LTD.
No.58 Jiyixiang Road, Langxia Street, Yuyao City,
Ningbo City, Zhejiang Province, China.
Tel: +86 18676591584

Importer:
Arzum Elektrikli Ev Aletleri San. ve Tic. A.Ş.
Defterdar Mah. Otakçılar Cad. Sinpaş Flatofis No: 78
İç Kapı No: 34 EYÜPSULTAN / İSTANBUL / TÜRKİYE
Made in China
Tel: +90 212 467 80 80
Fax: +90 212 467 80 00
www.arzum.com

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